

SIEMENS **Estimation of Contrast Agent Motion** from Cerebral Angiograms for Assessing Hemodynamic Alterations **Following Aneurysm Treatment by Flow Diversion**

T. Benz¹, M. Kowarschik^{1,2}, J. Endres³, P. Maday¹, T. Redel² and Nassir Navab¹

¹ Computer Aided Medical Procedures (CAMP), Technische Universität München (TUM), München, Germany ² Angiography & Interventional X-Ray Systems, Siemens AG, Healthcare Sector, Forchheim, Germany ³ Pattern Recognition Lab, Department of Computer Science, Friedrich-Alexander-Universität Erlangen-Nürnberg, Erlangen, Germany

Introduction

Flow diversion is an emerging treatment strategy for cerebral aneursyms that aims at restoring the healthy vessel anatomy by placing a tightly woven stent (flow diverter) across the opening of the aneurysm.





Once in place, the flow diverter (FD) redirects the blood flow away from the aneurysm and back into the parent vessel, thereby inducing flow stagnation and promoting aneurysm thrombosis.

Untreated aneurysm and gradual healing process after flow diverter treatment³

Objective: Angiographically decide whether the induced flow disruption is sufficient for aneurysm occlusion or whether additional therapeutic measures are needed.

Method





We recover dense contrast motion fields from pre/post pairs of Digital Subtraction Angiography (DSA) series by employing optical flow estimation based on conservation of mass [1] and second-order div*curl regularization* [2]. Metrics that capture hemodynamic changes induced by FD treatment are derived from the obtained motion fields in regions of interest (e.g., average speed inside the aneurysm before and after FD treatment).

Data Term	Regulari	ization Term
$I_t + \nabla I \cdot \vec{u} = 0 \qquad \qquad I_t + \nabla \cdot (I\vec{u}) = 0$	$\int \ \operatorname{div}(\vec{u})\ ^2 + \ \operatorname{curl}(\vec{u})\ ^2$	$\int \ \nabla \operatorname{div}(\vec{u})\ ^2 + \ \nabla \operatorname{curl}(\vec{u})\ ^2$
Constant brightness Conservation of mass	J_{Ω} First-order div-curl	J_{Ω} Second-order div-curl

Exemplary Results

Average Velocity Magnitude Pre vs. Post



Case A16: reduction in average speed inside aneurysm: **36%**



Reduction in average speed over whole cohort (n=17): 20%

On average (n=17 pairs of DSA series), we found a reduction in the average estimated speed inside the aneurysm of 20%.

References

[1] Wildes RP, Amabile MJ, Lanzillotto AM, Leu TS: Recovering estimates of fluid flow from image sequence data. CVIU 2000;80:246-266. [2] Corpetti T, Mémin É, Pérez P.: Dense estimation of fluid flows. PAMI 2002;24:365-380.

Figure Credits

¹© Sadasvian et al., Stroke, 40(3):952–958, 2009. ²© ev3 Endovascular, Inc., Plymouth, MN, U.S.A. ³© Mayo Clinic, Rochester, MN, U.S.A. All others courtesy of the authors.

Acknowledgements

We would like to give special thanks to Dr. Saruhan Çekirge, Hacettepe University Hospital, Ankara, Turkey, and Dr. Demetrius Lopes, Rush University Medical Center, Chicago, IL, U.S.A.

Contact

Tobias Benz, benzt@cs.tum.edu http://campar.in.tum.de



